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**Never send the entire DNA stock to the PGTB, only aliquots.**

### Sample preparation


⇒ Universal sequences must be added to your specific primers to enable hybridization of adapters and sequencing indexes in a 2nd PCR.

Amplicon PCR Forward Primer = 5' **ACACTCTTCCCTACACGACGCTCTCCGATCT**XXXXXXXXXXXXXXXXXXXX 3'  
5' **Universal sequence 1** – **Specific forward sequence 3'**

Amplicon PCR Reverse Primer = 5' **GTGACTGGAGTTCAGACGTGTGCTCTCCGATCT**XXXXXXXXXXXXXXXXXXXX 3'  
5' **Universal sequence 2** – **Specific reverse sequence 3'**

⇒ After PCR, check on an electrophoresis gel that it is correctly amplified to the expected sizes. Annotate your gel as shown here by specifying:

- type of size marker used
- the name of the samples or their position on the plate
- the expected fragment size

⇒ To deposit DNA in plates, use **96-well PCR plates that are transparent, skirted and non-break-away**. Each plate should contain **only one target**. If several targets are being analyzed in the project, separate them into different plates.

⇒ Identify each plate with **the acronym of the project declared in the project form** and a plate number: PROJECT X – PLATE 1.

⇒ Dispense (at minimal) 25 µL of the PCR products **in columns**. For PGTB controls, leave an empty well on the H line of each plate (plate 1 = H1; plate 2 = H2, etc.). If you have more than 12 plates, contact PGTB.

⇒ Heat seal the plates with suitable films (be careful not to deform them by overheating) or seal the plates with compatible cap strips. If you have used caps, wrap the plates in parafilm to keep the caps firmly in place.

### Depositing or sending samples

⇒ Complete the [sample submission form](#) to deposit samples at the PGTB or to send them by courier. Parcels are received from 9:00 a.m. to 5:30 p.m., Monday to Friday, at the following address:

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⇒ For shipping by courier, place ice packs at the bottom and sides of the package. Place the plates in the centre, separating them with cardboard to prevent them from being crushed. Fill any empty spaces with bubble wrap or cardboard so that the plates are securely held in place in the package.

⇒ PGTB will notify you by email when your samples have been received. **Any samples that do not meet the conditions specified in the email or in these requirements cannot be processed within the stated time frame. An additional charge may be applied or a new sample shipment may be requested.**

### Sample plate information

⇒ Once we have received the sample submission form, you will receive an email from PGTB with a link to directly access the file for submitting plate designs for your project on our RESANA platform.

⇒ Click on the link and the plate plan submission file will open in your web browser. Complete the plate plans following the instructions provided: **only letters/numbers/hyphens and underscores are permitted**.

⇒ You can copy/paste the plate designs using the keyboard shortcut 'Ctrl + V'.

⇒ **The file is automatically saved, so you do not need to do anything**. Once all the plate designs have been completed, close the file by clicking on the cross in the top right-hand corner.

	<b>Rédaction</b>	<b>Prescillia Alves Gomes - IE</b>	
	<b>Approbation</b>	<b>Zoé Delporte - AI</b>	