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	<b>Targeted metagenomics requirements using Illumina sequencer DNA SEND</b> (French version on page 1)	Instruction
		Version : 28.04.2026
		Pages : 2/2

**Never send the entire DNA stock to the PGTB, only aliquots.**

### Information on the conditions applied during locus-specific PCR

PGTB does not perform optimizations on PCR conditions: 35 cycles of PCR will be applied by default (unless otherwise requested by the user) and the hybridization temperature will be set at 55°C if no information has been communicated to us specifying it.

### DNA preparation

- ⇒ Quantify your DNA by **fluorescence**, then dilute it between 20 and 30 ng/μL. If the concentration is less than 10 ng/μL, inform the PGTB.
- ⇒ To deposit DNA in plates, use **96-well PCR plates that are transparent, skirted and non-break-away**.
- ⇒ Identify the plates using **the project acronym and the plate number**: PROJECT X - PLATE 1.
- ⇒ Dispense (at minimal) 30 μL of DNA **in columns**. For PGTB controls, **leave 2 empty wells on the H line of each plate (plate 1 = H1, H2; plate 2 = H2, H3, etc...)**. **Among the samples you are analysing, you must include a positive control to validate the PCR.**
- ⇒ Heat seal the plates with suitable films (be careful not to deform them by overheating) or seal the plates with compatible cap strips. If you have used caps, wrap the plates in parafilm to keep the caps firmly in place.

### Depositing or sending samples

- ⇒ Complete the [sample submission form](#) to deposit samples at the PGTB or to send them by courier. Parcels are received from 9:00 a.m. to 5:30 p.m., Monday to Friday, at the following address:

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- ⇒ For shipping by courier, place ice packs at the bottom and sides of the package. Place the plates in the centre, separating them with cardboard to prevent them from being crushed. Fill any empty spaces with bubble wrap or cardboard so that the plates are securely held in place in the package.
- ⇒ PGTB will notify you by email when your samples have been received. **Any samples that do not meet the conditions specified in the email or in these requirements cannot be processed within the stated time frame. An additional charge may be applied or a new sample shipment may be requested.**

### Sample plate information

- ⇒ Once we have received the sample submission form, you will receive an email from PGTB with a link to directly access the file for submitting plate designs for your project on our RESANA platform.
- ⇒ Click on the link and the plate design submission file will open in your web browser. Complete the plate plans following the instructions provided: **only letters/numbers/hyphens are permitted**.
- ⇒ You can copy/paste the plate designs using the keyboard shortcut 'Ctrl + V'.
- ⇒ **The file is saved automatically, so you do not need to do anything.** Once all the plate designs have been completed, close the file by clicking on the cross in the top right-hand corner.

	Rédaction	Prescillia Alves Gomes - IE	
	Approbation	Adline Delcamp – IE-RMQ	