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	Requirements for SSR and SNP genotyping (French version on page 1)	Instruction
		Version : 28.04.2026
		Pages : 2/2

Never send the entire DNA stock to the PGTB, only aliquots.

DNA preparation

- ⇒ For 'standard' DNA (extracted from tissue or recent material): quantify DNA by fluorescence, then dilute it to ~ 15 ng/μL. If concentrations are < 5 or > 25 ng/μL, the quality of the results may be compromised. Inform the PGTB.
- ⇒ For 'complex' DNA (degraded or ancient DNA, environmental DNA such as filtered water, faeces, etc.): do not dilute the DNA, proceed to the next step.
- ⇒ To deposit DNA in plates, use 96-well PCR plates that are transparent, skirted and non-break-away.
- ⇒ Identify each plate with the acronym of the project declared in the project form and a plate number: PROJECT X – PLATE 1.
- ⇒ Deposit 30 μL of DNA (diluted if necessary) into the plates, leaving one well empty without water or DNA on each plate as follows: plate 1 = H1 empty, plate 2 = H2 empty, ..., plate 12 = H12 empty, plate 13 = G1 empty, plate 14 = G2 empty, ..., plate 24 = empty G12, plate 25 = empty F1, etc.
- ⇒ Heat seal the plates with suitable films (be careful not to deform them by overheating) or seal the plates with compatible cap strips. If you have used caps, wrap the plates in parafilm to keep the caps firmly in place.

Depositing or sending samples

- ⇒ Complete [the sample submission form](#) to deposit samples at the PGTB or to send them by courier. Parcels are received from 9:00 a.m. to 5:30 p.m., Monday to Friday, at the following address:

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- ⇒ For shipping by courier, place ice packs at the bottom and sides of the package. Place the plates in the centre, separating them with cardboard to prevent them from being crushed. Fill any empty spaces with bubble wrap or cardboard so that the plates are securely held in place in the package.
- ⇒ PGTB will notify you by email when your samples have been received. Any samples that do not meet the conditions specified in the email or in these requirements cannot be processed within the stated time frame. An additional charge may be applied or a new sample shipment may be requested.

Sample plate information

- ⇒ Once we have received sample submission form, you will receive an email from PGTB with a link to directly access the file for submitting plate designs for your project on our RESANA platform.
- ⇒ Click on the link and the plate plan submission file will open in your web browser. Complete the plate plans following the instructions provided: only letters/numbers/hyphens are permitted.
- ⇒ You can copy/paste the plate designs using the keyboard shortcut 'Ctrl + V'.
- ⇒ The file is automatically saved, so you do not need to do anything. Once all the plate designs have been completed, close the file by clicking on the cross in the top right-hand corner.

	Rédaction	Adline Delcamp – IE-RMQ	
	Approbation	Erwan Guichoux - IR	